



# ACES™ ■ YebF Protein Export Kit

## AthenaES™ Complete Expression System

The ACES™ YebF Protein Export Kit uses the YebF protein to transport target proteins through the inner and outer membranes of *E. coli* cells to the extracellular medium in active form, allowing for the expression of toxic and difficult-to-fold proteins.



Athena Enzyme Systems™

888-892-8408

[www.athenaes.com](http://www.athenaes.com)

Application Manual

V 1.0



**ACES™**  
**AthenaES™ Complete Expression System**  
**YebF Protein Export Kit**

**Application Manual**  
**V. 1.0**

© Athena Environmental Sciences, Inc.

Turbo Broth™, Superior Broth™, Power Broth™, Hyper Broth™, Turbo Prime Broth™, Augmedium™, and ACES™ (AthenaES™ Complete Expression System™) are all trademarks of Athena Environmental Sciences, Inc. Athena Enzyme Systems™ is a division and trademark of Athena Environmental Sciences, Inc. AthenaES™ is a trademark of Athena Environmental Sciences, Inc. and the AthenaES™ logo is a registered logo of Athena Environmental Sciences, Inc.



**AthenaES®**

*a division of Athena Environmental Sciences, Inc.*

**Athena Enzyme Systems™**  
1450 South Rolling Road  
Baltimore, MD 21227

USA

T (MD): 410-455-6319  
T (USA): 888-892-8408  
F: 410-455-1155  
aesinfo@athenaes.com

# ACES™ YebF Protein Export Kit

## Application Manual

---

### **Table of Contents**

<b>Introduction</b>	<b>3</b>
<b>Principle of the Kit</b>	<b>4</b>
<b>Kit Components</b>	<b>5</b>
<b>Protocols:</b>	
<b>Preparation</b>	<b>6</b>
<b>Task 1: Subclone the Desired Gene Sequences into pAES40</b>	<b>8</b>
<b>Task 2: Express the YebF-Target Protein Fusion</b>	<b>8</b>
Basic Protocol for Task 2	8
Optimization Protocol for Task 2	9
<b>Task 3: Purification of the YebF Target Protein</b>	<b>11</b>
<b>Supplemental Protocols:</b>	
Rapid Transformation Protocol	12
Media Optimization Protocol	14
Determining Critical Factors for Expression Protocol	15
<b>References</b>	<b>16</b>
<b>Product Warranty</b>	<b>17</b>
<b>Product Limitations and License Information</b>	<b>17</b>
<b>Ordering Information</b>	<b>18</b>

## Introduction

*E. coli* is often the host strain of first choice for the production of recombinant proteins. Despite its long and successful history, several limitations remain when the desired protein can not be expressed in a functional state. This is often due to the inherent properties of expressing heterologous proteins in the cytoplasm. One means of overcoming these impediments is to express the protein such that it is secreted into an environment that is more conducive to correct folding and where proteolytic degradation is minimized, such as the periplasm or extracellular matrix. Periplasmic or extracellular protein export has been exploited for the production of a number of recombinant proteins.<sup>1</sup> The extracellular accumulation of a target protein provides several advantages. These include:

- **A Simplified Downstream Purification Scheme**

*E. coli* does not naturally export a significant number or amount of proteins to the extracellular matrix. Therefore, the level of host cell proteins as well as endotoxin and nucleic acids, contaminants that present considerable challenges to the purification of proteins destined for pharmaceutical use, would be significantly reduced.

- **Enhanced Biological Activity**

Export through the periplasmic space exposes proteins to a set of disulfide isomerases and foldases which facilitate correct protein folding as well as to an oxidizing environment which favors disulfide bridge formation.


- **Higher Product Stability and Solubility**

In addition to chaperones, there are fewer proteases in the periplasmic space and even fewer extracellular proteases.

However, *E. coli* does not normally secrete large numbers or quantities of proteins to the extracellular environment. Most secreted proteins reside in the periplasmic space. In order to export recombinant proteins to the culture medium a carrier protein is needed to effect transport across the outer membrane. One such protein shown to support protein export to the culture medium is YebF.<sup>2</sup> YebF has an unknown function, but is a *bona fide* extracellular protein. It also effectively transports both small and large prokaryotic and eukaryotic proteins to the extracellular medium in an active form. In this study, the *E. coli* expression vector pMS119 (Ap<sup>R</sup>, *ptac*),<sup>3</sup> was used to construct pYebFH6/MS. This plasmid expresses wild-type YebF protein under the control of an IPTG-inducible promoter, *ptac*, and with a C-terminal hexa-His affinity tag. Analysis of the subcellular localization of the YebFH6 protein after induction showed that the protein accumulated in the culture medium. To demonstrate that YebF could facilitate the export of other proteins, C-terminal fusions were made by inserting the coding sequences for mature alkaline phosphatase (*E. coli phoA*),  $\alpha$ -amylase (*Bacillus subtilis* X-23, *amy*) and human IL-2 between the C-terminal residue of YebF and the His tag. After induction all three proteins were found to accumulate in the culture medium, indicating that the YebF protein could effect extracellular transport of the fused protein. Importantly, cytoplasmic proteins did not leak into the medium. Therefore, YebF can be used to facilitate the extracellular export of recombinant proteins. Athena scientists have extended the work of Zhang *et al.*<sup>3</sup> by demonstrating that the YebF export function works in several strains of *E. coli* that are commonly used for the expression of heterologous proteins, including HB101, HMS174, BLR and TOP10.<sup>4</sup> The ACES™ YebF Protein Export Kit contains a user friendly vector, along with the accessory reagents needed for the successful production of a recombinant protein.

## Principle of the Kit

YebF is used as a carrier protein to direct the extracellular export of other recombinant proteins. This is achieved by fusing the desired protein to the C-terminus of YebF. After induction of expression, the fusion protein is exported to the culture medium where it accumulates. The fusion protein is recovered and the YebF portion removed using an endoprotease.



Blue - C-terminal amino acids of YebF

Orange - Enterokinase cleavage site

**Figure 1.** Plasmid map and MCS of pAES40

To use YebF as a carrier protein for the extracellular production of recombinant proteins, the gene encoding the desired protein is subcloned into the plasmid pAES40 (Fig. 1). The multiple cloning site for constructing the YebF-Target Protein fusion is shown in the vector drawing. The C-terminal two amino acids of YebF are the XhoI site with the remainder of the sequence written in the reading frame of YebF. An enterokinase proteolytic cleavage site (GAC GAT GAC GAT AAG) is between the multiple restriction sites used for subcloning and the end of YebF to permit removal of the YebF sequences after export. A hexa-His sequence is at the end of the MCS to provide an affinity tag if desired.

## Kit Components

ACES™ YebF Protein Export Kit Components		
Component	Amount	Catalog Number
pAES40	10 µg	0149-40
JM109 Competent Cells	2 x 200 µL	0151-JM109-C
Anti-YebF Antisera	0.5 mL	0313-1
LB Broth (Miller)	1 L Mix	0103
Turbo Broth™	1 L Mix	0104
Superior Broth™	1 L Mix	0105
Power Broth™	1 L Mix	0106
Hyper Broth™	1 L Mix	0107
Glucose M9Y	1 L Mix	0108
Glucose Nutrient Mix	21 g	0109
Augmedium™	25 mL Mix	0123
Inducer Solution A	1 mL	0152
Inducer Solution B	100 mL	0153
Secretion Enhancer Solution A	60 mL	0154
Secretion Enhancer Solution B	60 mL	0155
Reagents needed but not provided: Glycerol, LB plates supplemented with 100 µg/ml ampicillin		

Solution Information	
Component	Reagent
Inducer Solution A	1 M IPTG
Inducer Solution B	20% Lactose
Secretion Enhancer Solution A	10% Glycine
Secretion Enhancer Solution B	10% Triton X-100

## Protocols

### Preparation: Media

1. Dissolve the contents of each of the media packets in deionized water as directed on the individual packets.
2. Add 4 mL of glycerol to the Turbo Broth™ and Power Broth™ solutions.
3. Dispense desired volume into appropriate bottles or flasks. (We recommend 2 x 500 mL glass bottles.)
4. Autoclave at 121°C for 20 min. The autoclaved media without antibiotics are stable for at least 6 months at 4°C.
5. Dissolve the contents of the Glucose Nutrient Mix in 100 mL deionized water and filter sterilize using a 0.2 µm filter.
6. Add 50 mL of the sterile Glucose Nutrient Mix to 1 liter of Hyper Broth™ and 20 mL to 1 liter of Glucose M9Y using aseptic technique.
7. Add sterile antibiotics as needed.

### Preparation: Medium Supplements

#### 1. Augmedium

- 1.1. Dissolve the contents of the Augmedium™ container in 25 mL deionized water.
- 1.2. Filter sterilize using a 0.2 µm filter.
- 1.3. Aliquot 5 mL portions into sterile 15 mL conical tubes and store at -20°C. The stock will be a 50x solution.
- 1.4. For use, snap thaw at 37°C and store at 4°C for no more than five (5) days.

#### 2. Inducer Solution A (1M IPTG)

- 2.1. Dilute to 10 mM by adding 5 µl to 495 µl sterile water.
- 2.2. Store at -20°C.

#### 3. Inducer Solution B (20% Lactose)

- 3.1. Make 20% stock solution by dissolving contents of packet into 100 mL of deionized water.

3.1. Filter sterilize using 0.2 µm filter.

3.2. Store at 4°C.

#### 4. Secretion Enhancer Solution A (10% Glycine)

4.1. Make 10% stock solution by dissolving contents of packet into 60 mL deionized water.

4.2. Filter sterilize using 0.2 µm filter.

4.3. Store at 4°C.

#### 5. Secretion Enhancer Solution B (10% Triton X-100)

5.1. Use as is to a final concentration of 1%.

5.2. Store at 4°C.



**Task 1: Subclone the Desired Gene Sequences into pAES40**

Use any of the restriction sites located downstream of the enterokinase cleavage site to insert the coding sequence for the target protein. Be sure to design the subcloning such that the reading frame of the target protein matches that of YebF.

**Task 2: Express the YebF-Target Protein Fusion**

The protocol presented below is the basic method for producing and recovering proteins exported using the YebF Protein Export System. Suggested alternative methods to this basic scheme are noted. Briefly, a Lac operon-based induction protocol is employed. Inducer Solution A or B can be used as the inducer. Inducer Solution A is used initially at 50  $\mu$ M and can be increased to 1 mM depending on the viability of the host strain and production levels of the target protein. Inducer Solution B is used at 2%. For the YebF Protein Export System, induction times should be greater than 16 hours to allow the exported protein to accumulate in the culture medium. Any rich medium can be employed.

**• Basic Protocol for Task 2: (For 1 Liter, Adjust Volumes as Needed)**

1. Introduce the plasmid construct from Step 1 into the desired production strain.
  - 1.1. Supplemental Protocol 1 gives a protocol for the rapid transformation of any bacterial vector into any *E. coli* strain.
  - 1.2. Streak purify two to eight transformants and verify the integrity of the plasmid construct.
  - 1.3. Store the strains accordingly.
2. Use a single colony to inoculate a 50 ml starter culture supplemented with 100  $\mu$ g/ml ampicillin.
  - 2.1. Turbo Broth™ or Turbo Prime Broth™ (AthenaES™ Cat. Nos. 0104 and 0110, respectively) are recommended but any rich medium will suffice.
  - 2.2. Incubate at 37°C overnight.
3. Inoculate 1 liter of medium supplemented with 100  $\mu$ g/ml ampicillin with the overnight culture and incubate at 30°C until the OD<sub>600</sub> reaches 0.8-1.0.
  - 3.1. The use of baffle bottomed flasks is strongly recommended to provide sufficient aeration.
  - 3.2. Turbo Broth™ or Turbo Prime Broth™ are recommended but any rich medium will suffice.

4. Induce expression by adding 50 µl Inducer Solution A (IPTG ) (to 50 µM), or 100 ml Inducer Solution B (Lactose) (to 2%).
  - 4.1. Note: Inducer Solution B can only be used if the host strain is wild-type for lacZ and lacY.
5. Incubate at 30°C for 20-24 hours.
6. Remove the cells by centrifugation at 5,000 xg for 30 min.
7. Add 472 g ammonium sulfate to the culture supernatant and dissolve thoroughly.
8. Chill on ice for 2 hours with gentle stirring.
9. Collect the precipitated protein by centrifugation at 15,000 xg for 30 min. at 4°C.
10. Suspend the pellet in 25 ml of the desired buffer, filter with 0.45 µm cartridge or depth filter, and proceed with subsequent purification steps as needed.

- **Optimization Protocol for Task 2:**

1. Introduce the plasmid construct from Step 1 into each of three to six host strains.
  - 1.1. The most commonly used strains for the production of recombinant proteins can be used. See Supplemental Protocol 1 for a simple rapid transformation protocol. It is recommended that the parent pAES40 plasmid be introduced into each strain to serve as a control in expression experiments.
2. Select the strains with the highest level of export of the fusion protein.
  - 2.1. Inoculate 5 ml of Turbo Broth™ (AthenaES™ Cat. No. 0104) supplemented with 100 µg/ml ampicillin with a single colony of each strain harboring the pAES40 derivative. Incubate the culture at 37°C overnight.
  - 2.2. Use the overnight culture to inoculate duplicate cultures of 25 ml fresh Turbo Broth™ supplemented with 100 µg/ml ampicillin. Incubate the cultures at 30°C until the absorbance at OD<sub>600</sub> reaches 0.8. Remove a 1 ml sample, pellet the cells and store the culture supernatant at 4°C. (Note: The cell pellet can be reserved by storing at -20°C.)
  - 2.3. Induce expression by adding 125 µl diluted Inducer Solution A (IPTG ) to (50 µM), or 2.5 ml Inducer Solution B (Lactose) (to 2%).
    - 2.3.1. Dilute Inducer Solution A (1 M) by adding 5 µl to 495 µl sterile water (final concentration 10 mM). Add 125 µl of the diluted solution to each culture. Store the remainder of the diluted solution at -20°C for use in later experiments.

Critical Factors that Affect Accumulation of the Target Protein						
Culture	IPTG (mM)	Lactose (%)	Glycine (%)	Triton X-100 (%)	Augmedium™ (x)	Temp.
1	1	0	1	0	0	37°C
2	0.05	2	1	0	0	27°C
3	1	2	1	1	1	37°C
4	1	0	0	0	1	27°C
5	1	2	0	0	0	37°C
6	1	2	1	0	1	27°C
7	0.05	0	0	0	0	27°C
8	1	0	0	1	1	37°C
9	0.05	2	1	1	0	37°C
10	1	0	1	1	0	27°C
11	0.05	2	0	1	1	27°C
12	1	2	0	1	0	27°C
13	0.05	0	0	1	0	37°C
14	0.05	0	1	0	1	37°C
15	0.05	0	1	1	1	27°C
16	0.05	2	0	0	1	37°C

**Table 1.** Experimental design matrix for determining critical factors that affect accumulation of a target protein.

- 2.3.2. To the duplicate culture, add 2.5 ml of Inducer Solution B. (Note: In some cases expression may be enhanced by using 2% lactose in place of IPTG as the inducer. This is suitable only for strains which are wild-type for lacZY.)
- 2.4. At 3-6 hours post-induction remove a 1 ml sample and process as in step 2.2.
- 2.5. At 20-24 hours post-induction remove a 1 ml sample and process as in step 2.2.
- 2.6. Remove the cells by centrifugation at 5,000 xg for 20 min and store the culture supernatant at 4°C.
- 2.7. Analyze the culture supernatants by immunoblot, ELISA or functional assay. Do not use a stained gel to visualize the expression of the YebF fusion as the accumulation levels in a non-optimized system may be below the detection limits of SDS-PAGE. The immunoassays can be done with anti-Target Protein, anti-YebF (AthenaES™ Cat. No. 0313), or anti-His tag antibodies.
- 2.8. Select the strain that produced the most YebF-Target Protein per ml of culture.
3. Determine the optimum medium for production.

- 3.1. Perform a media screen as described in Supplemental Protocol 2 and select the medium that yields the highest level of Target Protein per ml.
4. Optimize the culture conditions for production of the fusion protein.
  - 4.1. Determining the optimal culture conditions is an iterative process. First, a two-level fractional factorial experimental design is used to identify the critical factors affecting the accumulation of the Target Protein.<sup>5</sup> Table 1 lists the conditions to be tested. Use an immunoassay or functional assay to quantify the level of expression and select the conditions that yield the highest level of production. Supplemental Protocol 3 provides a method for interpreting the data. Further optimization is then done by modifying those parameters identified as critical to accumulation of the Target Protein.

### Task 3: Purification of the YebF-Target Protein

Using the strain and culture conditions identified in Task 2, perform a production run using 0.5 to 1 liter sized cultures. It is recommended that baffle bottomed flasks be used to ensure adequate aeration. The YebF-Target Protein is recovered from the culture medium as follows:

1. After 20-24 hours post-induction, remove the cells by centrifugation at 5,000 xg for 20 min.
2. For each liter of medium, add 472 g of ammonium sulfate.
3. Allow the salt to completely dissolve and then chill the medium on ice for 2 hours with gentle stirring.
4. Dispense the mixture into centrifuge bottles and collect the precipitated protein by centrifugation at 15,000 xg for 30 min at 4°C.
5. Decant the liquid and suspend the pellets in 25 ml of a buffer suitable for the next purification step.
  - 5.1. To ensure complete removal of any residual cells or debris, filter the solution using a 0.45 µm cartridge or depth filter.
6. Dialyze against the enterokinase (EK) buffer to ensure removal of the ammonium sulfate and process as needed.
7. Cleave YebF from the Target Protein using EK.
  - 7.1. Note: The suggested starting concentration for EK is 1 Unit of enzyme/25µg of fusion protein. More or less enzyme may be required for any given Target Protein.
8. Purify the Target Protein accordingly.

## Supplemental Protocols

### Supplemental Protocol 1: Rapid Transformation Protocol (Kit available: AthenaES™ No. 0156)

#### 1. Materials:

- 1.1. 2x TSS (DO NOT store or use near flame, DMSO is flammable).
- 1.2. *E. coli* strains to be transformed.
- 1.3. Tryptic Soy Agar (TSA) plates or comparable non-antibiotic containing plate medium.
- 1.4. Sterile microcentrifuge tubes, 1.5mL (or other convenient tube).
- 1.5. LB medium.
- 1.6. SOC medium.
- 1.7. Antibiotic-containing plate medium.

#### 2. Methods:

- 2.1. Streak strain(s) on non-antibiotic plate medium and incubate at 37°C overnight. (This protocol MUST be performed with fresh, overnight bacterial colonies.)
- 2.2. Dispense 0.1 mL LB medium into microfuge tubes using sterile technique.
- 2.3. Pick 4 colonies using a 1µL loop or sterile toothpick and resuspend in the 0.1 mL LB using sterile technique.
- 2.4. Add 0.1 mL 2x TSS (DO NOT FLAME) and mix well.
- 2.5. Incubate on ice for 15 minutes. Once chilled, do not allow the cells to warm above 14°C.
- 2.6. Add 100 ng plasmid DNA and incubate on ice for 20 min.
- 2.7. Heat shock in a 42°C water bath for 1 min.
- 2.8. Add 1 mL SOC medium using sterile technique.
- 2.9. Incubate at 37°C for 30 min.
- 2.10. Plate 0.1 mL of transformation mix on antibiotic-containing plate medium and incubate overnight at 37°C for 30 min.

- 2.11. Streak purify 2 or 3 colonies on appropriate antibiotic plates.
- 2.12. Prepare master cell bank(s) or cryostock(s) of new strain(s).

## Supplemental Protocol 2: Media Optimization Protocol

### 1. Materials

- 1.1. 25 mL of each culture medium in 250 mL baffle bottomed flasks.
- 1.2. Reagents for Immunoblot or Functional Assay

### 2. Methods

- 2.1. Inoculate a single colony of the recombinant strain into 10 ml of LB Broth in a shake flask with a baffle bottom. Incubate at 37°C overnight.
- 2.2. Inoculate 25 ml of each of the six media with 2 ml of the overnight culture. Incubate the cultures at 30°C until the OD<sub>600</sub> reaches 0.8.
- 2.3. Remove a 1 ml sample (“pre-induction”), harvest the cells and store the culture medium at 4°C for later analysis.
- 2.4. Add the inducer and continue incubating for 22 hours.
  - 2.4.1. Add 125 µl of 1:100 diluted Inducer Solution A (final IPTG concentration of 50 µM) or 2.5 ml of Inducer Solution B (final lactose concentration of 2%).
- 2.5. Remove a 1 ml sample (“post-induction”) and process as in step 2.3.
- 2.6. Remove the cells from the remaining culture medium by centrifugation at 5,000 xg for 20 min and store the supernatant at 4°C.
- 2.7. Determine the relative level of Target Protein in the pre-induction and post-induction samples for each culture condition using an immunoblot, (detection with anti-YebF, anti-His Tag or anti-Target Protein antibodies), or a functional assay.
- 2.8. Select the medium yielding the highest level of Target Protein accumulation for further work.

**Supplemental Protocol 3: Determining the Critical Factors for Expression.**

1. Prepare a spreadsheet with 16 rows corresponding to culture conditions 1 through 16 and 6 columns corresponding to each of the factors tested and the solutions as shown in the figure below.
2. Enter the value (i.e., enzyme activity, mass, etc.) obtained for each culture condition into the respective cell in the row. For any given condition each factor will have the same value entered.
  - 2.1. Note: Numeric descriptors for qualitative assessments will also work, but with less accuracy.)
3. Calculate the sum of protein produced for each factor when the factor was present in the solution.  $\text{Sum}_{\text{Present}}$ 
  - 3.1. Note: For IPTG consider the 1 mM level “Present.”
4. Calculate the sum of protein recovered for each factor when the factor was absent from the solution.  $\text{Sum}_{\text{Absent}}$ 
  - 4.1. Note: For IPTG consider the 0.05 mM level “Absent.”
5. Calculate the difference between the Present and Absent and divide by 8 for each factor.  $\text{Relative Effect} = \text{Sum}_{\text{Present}} - \text{Sum}_{\text{Absent}} / 8.$
6. Compare the Relative Effect numbers obtained.
  - 6.1. A positive number indicates a positive effect on accumulation.
  - 6.2. A negative number indicates no effect on accumulation.
  - 6.3. The larger the positive number the greater the effect of the given factor.

Critical Factors that Affect Accumulation of the Target Protein						
Culture	IPTG (mM)	Lactose (%)	Glycine (%)	Triton X-100 (%)	Augmedium™ (x)	Temp. (°C)
1						
2						
.						
.						
.						
15						
16						
$\text{Sum}_{\text{Present}}$						
$\text{Sum}_{\text{Absent}}$						
Rel. Effect						

**Table 2.** Example analysis table used to determine the factors critical to protein accumulation.



## References

1. Choi, J. H., and Lee, S. Y. 2004. Secretory and extracellular production of recombinant protein using *Escherichia coli*. *Appl. Microbiol. Biotechnol.* 64:625-635.
2. Zhang, G., Brokx, S. and Weiner, J. H. 2006. Extracellular accumulation of recombinant protein fused to the carrier protein YebF in *Escherichia coli*. *Nat. Biotech.* 24:100-104.
3. Strack, B., Lessel, M., Calendar, R. and Lanka, E. 1992. A common sequence motif, -E-G-Y-A-T-A-, identified within the primase domains of plasmid-encoded I- and P-type DNA primases and the  $\alpha$  protein of the *Escherichia coli* satellite phage P4. *J. Biol. Chem.* 267:13062-13072.
4. Broedel, Jr., S. E. and Papciak, S. M. 2007. ACES™ Signal Sequence and YebF Expression Systems. Athena Environmental Sciences, Inc., Technical Brief, December 2007, [http://athenaes.com/osc/TechBrief\\_ACESSignalSeq\\_Web.pdf](http://athenaes.com/osc/TechBrief_ACESSignalSeq_Web.pdf).
5. Montgomery, D. C. 2001. *Design and Analysis of Experiments*. John Wiley & Sons, Hoboken, NJ, ISBN No. 0-471-31649-0.

## **Technical Assistance**

The scientific staff of the Athena Enzyme Systems™ are specialists in the expression of recombinant proteins in microbial systems. They have extensive expertise in all aspects of protein expression from the construction of expression vectors to the commercial production of recombinant proteins. No matter what your question, please feel free to ask them for help. A technical support scientist can be reached at [support@athenaes.com](mailto:support@athenaes.com).

## **Product Use Limitations**

The ACES™ YebF Protein Export Kit and Vector were designed and are sold for research use only. None of the kit components should be used for human diagnosis or drug use or administered to humans unless expressly cleared for that purpose by the appropriate regulatory authorities in the country of use. All due care and attention should be exercised in the handling of the materials contained in the kit.

The plasmids and other components contained in the ACES™ YebF Protein Export Kit are covered under international patents and patent applications exclusively licensed by the University of Alberta to Athena Environmental Sciences, Inc. The use of these reagents requires a license.

- Academic and Not-for-Profit entities are required to execute a non-commercial use license. This agreement is a non-fee bearing license which grants the user the rights to use the plasmid and other kit components for research purposes only and restricts the user from disseminating the plasmids to other researchers without the expressed written consent of AthenaES.
- Commercial users are required to execute a commercial evaluation license agreement. The agreement is a non-fee bearing license which grants to the user the right to use the plasmid and other kit components for research purposes only for the period of one (1) year after which a commercial use license is required.

Copies of both types of license agreements are available at [www.athenaes.com](http://www.athenaes.com).

## **Product Warranty**

AthenaES™ guarantees the quality and performance of the media and reagents contained in this kit for the cultivation of *E. coli*. The suitability of a medium formulation or additive for a particular use is the responsibility of the end user. No guarantee is made that a given protein will be expressed when applying this kit. AthenaES™ will replace the product free of charge if it does not conform to the stated specifications. Notice for replacement must be received within 60 days of opening the product.

## Ordering Information

### To place an order:

Phone: 1-888-892-8408 Email: [media@athenaes.com](mailto:media@athenaes.com)  
 Fax: 410-455-1155 Website: [www.athenaes.com](http://www.athenaes.com)

Or visit our website to order through one of our international distributors.

### When placing an order, please provide the following:

- Institution name and customer service account
- Purchase order number
- Catalog number(s) or names of products and quantity of item(s)
- Billing and shipping address
- Contact name and telephone number

### Delivery:

Telephone orders received Monday through Friday before 12 noon will be shipped that day. All other orders will be shipped the next business day, unless otherwise stipulated.

#### ACES™ Accessory Products Information

Catalog Number	Product	Size
0149-25	pAES25	10µg
0149-30	pAES30	10µg
0149-31	pAES31	10µg
0149-32	pAES32	10µg
0149-33	pAES33	10µg
0149-34	pAES34	10µg
0149-35	pAES35	10µg
0149-40	pAES40	10µg
0150-1	Primer A	250pmoles
0150-2	Primer B	250pmoles
0151-JM109-C	JM109 Competent Cells	2 x 200µL
0151-	Strain Stab	Stab
0152-1	Inducer Solution A	1mL
0152-5	Inducer Solution A	5 x 1mL
0153	Inducer Solution B	500mL
0154	Secretion Enhancer A	500mL
0155	Secretion Enhancer B	500mL
0156	Rapid Transformation Kit	1 kit
0157	2x TSS	5 x 1 mL
0313-1	Anti-YebF Antisera	0.5mL