

Protein Production and Purification




“Why the medium composition counts.”

Simple Solutions
for Complex Proteins



The *E. coli* System



Advantages:

- Well developed science and tools with a 40 year history of successes
- Inexpensive
 - Short development timelines
 - Inexpensive tools
 - Low cost production systems


Disadvantages:

- Endotoxin contamination
- Heterologous proteins often accumulate as insoluble products in cytoplasm – “Inclusion bodies”
- Export to periplasm is possible but rare for secretion to the medium
- No post-translational glycosylation

Simple Solutions
for Complex Proteins



Variables in Protein Production



Adapted from Gustafsson et al 2012. Protein Expression and Purification. 83:37-46.

Simple Solutions
for Complex Proteins



Variables in Protein Purification

Extraction

Intra- or Extracellular
Host contaminants
Process contaminants

Purification

Physical-chemical properties
Specifications
Protein stability
Protein-protein interactions
Means of measuring
Host contaminants

Stabilization

Degradation pathways
Trace inhibitors

- Empirical science
- Highly interconnected
- Choice of host-vector system
 - Type of extraction process
 - Post-translational modifications

Simple Solutions
for Complex Proteins



E. coli Expression Vectors

- Inducible, cytoplasmic accumulation
 - Lac Operon: promoters, T7, T5, *tac*, *trc*, *lacP*
 - Arabinose Operon
 - Tryptophan Operon
 - Lambda Repressor
- With Affinity Tags
 - His, MalE, GST, Trx, Flag, and others
- Periplasmic Expression
 - Signal Sequences: OmpA, PelB, SpA, PhoA, and others

The History of LB Broths

- Originally developed in the 1950's to cultivate *E. coli*.
- Comes in three variations
 - Miller, Lennox and Luria
 - Differ in NaCl content: 10, 5 and 0.5 g/L, respectively
- Period in time when optimum growth conditions were not known.
- Long before recombinant proteins were produced in *E. coli*.
- And they worked for the physiological and genetic experiments of the time.
- *However,*

The LB Broths

- Have no added carbon source.
- Are not buffered.
- No added phosphate, sulfate, or potassium

Composition of LB Broths	
Yeast Extract	5 g/L
Casein Hydrolysate	10 g/L
NaCl	10 g/L (Miller) 5 g/L (Lennox) 0.5 g/L (Luria)

Not designed or intended for production of recombinant proteins.

Better Media Formulations

Biomass yield of *E. coli* grown in six different medium.

Medium	Biomass Yield (g/L)
LB (Miller) Broth	10
Glucose M9Y	16
Hyper Broth™	36
Power Broth™	24
Superior Broth™	18
Turbo Broth™	30

E. coli strain JM109 was grown in 100 ml shake-flask cultures in each medium at 37°C for 16 h.

- Not all proteins express well in *E. coli*.
- 20 years of helping clients overcome this limitation.
- Discovered that the carbon and nitrogen source can make a big difference.

Simple Solutions
for Complex Proteins



Medium Composition *Makes a Difference*





Figure. SDS-PAGE analysis of total protein from each strain in Table 1. Samples were prepared as described in the Materials and Methods section. Panel A - reference strain without a recombinant protein; Panel B to G - MalE, GST, GFP, I278, TesA, LypA, respectively. Arrows denote the location of the respective protein. Marker proteins are shown to the left and right of each set of cellular proteins. From left to right in each panel are samples from cells grown in LB (Miller), Glucose M9Y, Hyper Broth™, Power Broth™, Superior Broth™ and Turbo Broth™.

Simple Solutions
for Complex Proteins



Medium Composition *Makes a Difference*

**Mammalian protein produced in a
400 liter fermentor**





Lane 1 = MW marker
Lane 2 = Un-induced fermentor
Lane 3 = Fermentor harvest, IPTG
Lane 4 = UGA's Auto-induction
Lane 5 = Hyper Broth auto-induction
Lane 6 = Power Prime auto-induction
Lane 7 = Overnight Express™ auto-ind.
Lane 8 = MW marker.

Overnight Express™ Autoinduction Medium is a trademark of EMDMillipore.

Simple Solutions
for Complex Proteins



Not Just More Protein, But More Active Protein



Simple Solutions
for Complex Proteins



Recommendations

- Richer media formulations will yield more recombinant protein than LB Broths.
- The best medium to produce a protein in a given host-vector system should be empirically determined.
- A simple screen will make all the difference.



Simple Solutions
for Complex Proteins



Where to Buy

- Screening Kits are available at Athena or one of our distributors.
- Visit www.athenaes.com/Expression.php
 - Look under “Protein Expression Media”



Media Optimization Kit

Available in standard and animal product free versions.

Simple Solutions
for Complex Proteins



References

1. Strobel, R. J. and Sullivan, G. R. 1999. Experimental design for improvement of fermentations. *In* Manual of Industrial Microbiology and Biotechnology, 2nd Ed., eds., Demain, A. L. and Davies, J. E. ASM Press, Washington, DC, pg 80-93.
2. Lee, S. Y. 1996. High Cell-density culture of *Escherichia coli*. *Tibtech* 14:98-105.
3. Kleman, G. L. and W. R. Strohl. 1994. Development in high cell density and high productivity microbial fermentation. *Current Opinion in Biotech.* 5:180-186.
4. Broedel, S. E., Papciak, S. M. and Jones, W. R. 2001. The Selection of Optimum Media Formulations for Improved Expression of Recombinant Proteins In *E. coli*. Athena Technical Bulletin, http://www.athenaes.com/tech_brief_optimum_media.php
5. Luria, S. E. and J. W. Burrous. 1955. Hybridization between *Escherichia coli* and *Shigella*. *J. Bacteriol.* 74:461-476.
6. Lennox, E. S. 1955. Transduction of linked genetic characters of the host by bacteriophage P1. *Virology.* 1:190-206.
7. Luria, S. E., J. N. Adams, and R. C. ting. 1960. Transduction of lactose-utilizing ability among strains of *E. coli* and *S. dysenteriae* and the properties of the transducing phage particles. *Virology.* 12:348-390.
8. Miller, 1972. Experiments in Molecular Genetics. Cold Spring Harbor Laboratory. Cold Spring Harbor, NY.